



Power-PRO TL

powered ambulance cot

stryker®

Reduce the risk of injuries when raising and lowering



Stryker SA
Cit -Centre
Grand-Rue 90
1820 Montreux
Switzerland

t: +41 21 966 12 01
f: +41 21 966 12 00

www.stryker.eu

Standard Features

- Manual back-up release handle
- Automatic high-speed retract
- Battery-powered hydraulic lift system
- Preventative maintenance hour meter
- Battery pack power level indicator
- Colour-coded controls
- High visibility powder-coated frame
- Lightweight, rugged aluminium construction
- Integrated bumper system
- Perforated litter surface
- One-hand release, fold down side rails
- One-hand release, infinite positioning, pneumatically assisted backrest
- Oversized wheels with sealed caster and wheel bearings
- Sealed bolster mattress
- Two lap belts and one four-point shoulder restraint
- Dual wheel locks
- Power washable
- Shock or flat leg positioning
- Head end steer-lock system
- Retractable head and foot sections

Optional Features

- Heavy duty two- and three-stage IV poles (patient right or left)
- Restraint extender (91.4 cm/36 in.)
- Corner poles (removable)
- Removable O₂ bottle holder
- Base storage net
- Pocketed backrest storage pouch
- Equipment hook
- Four-corner wheel locks
- SMRT Power System (12V DC, 120V AC and 240V AC options available)
- SMRT Charger mounting bracket
- Knee gatch
- Defibrillator platform
- Transfer Flat
- Pedi-Mate restraint package

Specifications

Model Number	6550
Height¹ (infinite height positioning between lowest and highest position)	
Highest Position	96.5 cm (38 in.)
Lowest Position	49.5 cm (19.5 in.)
Length	
Standard	194.3 cm (76.5 in.)
Minimum	139.7 cm (55 in.)
Width	58 cm (23 in.)
Weight²	66 kg (145 lb.)
Wheels	
Diameter	15 cm (6 in.)
Width	5 cm (2 in.)
Articulation	
Backrest	0 – 75°
Shock Position	+15°
Optional Knee Gatch	30°
Maximum Weight Capacity	318 kg (700 lb.)
Minimum Operators Required	
Occupied Cot	2
Unoccupied Cot	1
Recommended Fastener System	
Floor Mount ³	Model 6385, 6386, 6387

Warranty

- One-year parts, labour and travel or two-years parts
- Lifetime on all welds*
- Seller's warranty hereunder shall not apply if: (i) Goods are not used in accordance with instructions or are used for a purpose not indicated on the labelling, (ii) any repairs, alterations or other work has been performed by Buyer or others on such item, other than work performed with Seller's authorization and according to its approved procedures; or (iii) the alleged defect is a result of abuse, misuse, improper maintenance, accident or the negligence of any party other than Seller. The warranty set forth herein is conditioned upon proper storage, installation, use and maintenance in accordance with applicable written recommendations of Seller. The warranty furnished hereunder does not extend to damage to Goods purchased hereunder resulting in whole or in part from the use of components, accessories, parts or supplies not furnished by Seller.

Seller's sole obligation shall be to repair or replace, at Seller's option, any defective component or item and pay transportation expenses for such replacement.

*7-year service life.

¹ Height measured from bottom of mattress, at seat section, to ground level.

² Cot is weighed with one battery pack, without mattress and restraints.

³ Designed to be compatible with competitive ambulance cots. Floor plate and floor track mounting configurations available.

The Power-PRO TL is designed to be compatible with competitive cot fastener systems.

Stryker reserves the right to change specifications without notice.

This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Stryker**. All other trademarks are trademarks of their respective owners or holders.

The products listed above are CE marked according to the Medical Device Directive 93/42/EEC

Electrical and Mechanical Safety Certifications



Literature Number: MTX14901EN Rev. D
MTX-RRD/GS 01/11

Copyright © 2011 Stryker



MTX 14901EN