

stryker®



XPS

expandable patient surface

XPS

expandable patient surface



The all-in-one solution for a variety of patients and environments.

XPS is there when you need it.

XPS provides an expanded patient surface area that can easily be retrofitted with compatible cots. XPS is adjustable with 7 locking positions and includes a wider mattress that reduces transfer gap and is designed with patient comfort in mind. This solution helps address growing obesity trends and supports a variety of patients and environments.

Cots with XPS meet current tip stability standards for both patient and provider safety¹. Engineered for durability and cleanability, with a 6061-T6 aluminium core and over-mold design.

Standard Features:

- Increased patient surface area
- Adjustability for patients and environments
- 7 locking positions
- Integrated into cot – Always there
- Easily retrofitted to existing compatible cots
- Enhanced safety built into release handles
- Durable aluminium over-mold design
- Mattress design reduces transfer gap
- Enhanced patient comfort
- Compliant with tip stability and dynamic crash certifications¹




Release Handles

Enhanced safety built into release handles.

¹Certified to IEC 60601-1 for Power-PRO XT and Power-PRO TL.
BS-EN 1789 for Power-PRO XT, Power-PRO TL, and Performance-PRO XT.


7 Locking Positions

Between 10° and 52°



Increased Surface Area

From 58 cm to 84 cm



Option or Retrofit

XPS Compatible Cots

Power-PRO XT Model 6500/6506

Power-PRO TL Model 6550

Performance-PRO XT Model 6085/6086



Safety

Compliant with tip stability and dynamic crash certifications.

Durability

Engineered for durability with an aluminium core.

Always There

Integrated into cot for accessibility on every call.

Reduced Transfer Gap

Mattress design reduces transfer gap.

XPS Specifications

Includes

2 XPS Siderails
XPS Mattress

Weight¹ 8 lb (4 kg)

7 Locking Positions Between 10°-52°

Surface Area Expansion

Locking Angle	Width - Total Surface Area
10°	23 in (58 cm)
17°	25 in (64 cm)
24°	27 in (69 cm)
31°	29 in (74 cm)
38°	30 in (76 cm)
45°	32 in (81 cm)
52°	33 in (84 cm)
Height	10 in (25 cm)
Length	30 in (76 cm)
XPS Mattress Width²	23 in (58 cm)
XPS Compatible Cots³	
Power-PRO XT	Model 6500/6506
Power-PRO TL	Model 6550
Performance-PRO XT	Model 6085/6086

¹Additional weight compared to standard siderails not including mattress

²Width is measured at widest point (Standard bolster mattress 19 in / 48 cm)

³Certified to IEC 60601-1 for Power-PRO XT and Power-PRO TL.

BS-EN 1789 for Power-PRO XT, Power-PRO TL, and Performance-PRO XT.

Patents pending.

Stryker reserves the right to change specifications without notice.

All numbers rounded to nearest whole value.

This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Performance-PRO XT**, **Power-PRO TL**, **Stryker** and **XPS**. All other trademarks are trademarks of their respective owners or holders.

The products listed above are CE marked according to European Medical Device Directive.

This material is not intended for distribution outside the EU and EFTA areas.



Literature number: MTX81319EN Rev B.3.

MTX-6987-13-RRD/GS 2013/06

Copyright © 2013 Stryker

Reconstructive

Hips
Knees
Trauma & Extremities
Joint Preservation
Orthobiologics

Medical & Surgical

Power Tools & Surgical Accessories
Image Guided Navigation
Endoscopy & Arthroscopy
Integrated Communications
Beds, Stretchers & EMS
Sustainability Solutions

Neurotechnology & Spine

Craniofacial
Interventional Spine
Neurosurgical, Spine & ENT
Neurovascular
Spinal Implants



Stryker Medical

3800 E. Centre Ave.
Portage, MI 49002 U.S.A.
t: 269 329 2100
f: 269 329 2315
toll free: 800 784 4336
www.ems.stryker.com



Stryker France SAS

ZAC - Avenue de Satolas Green
Pusignan, 69330
France
t: +33 4 72 45 36 00
f: +33 4 72 45 36 99
www.stryker.fr

Distributed by:

Stryker SAS
Cité Centre
Grand-Rue 90
1820 Montreux
Switzerland
t: +41 21 966 12 01
f: +41 21 966 12 00