

Reduce the risk of injuries when on the stairs



Stryker SA
Cité-Centre
Grand-Rue 90
1820 Montreux
Switzerland

t: +41 21 966 12 01
f: +41 21 966 12 00

www.stryker.eu

Standard Features

- High visibility powder-coated frame
- Colour-coded controls
- Extendable foot end lift handles
- Locking rear lift handles
- Molded hand grips
- Lightweight, rugged aluminium construction
- Oversized rear wheels with sealed bearings
- Dual wheel locks
- Compact storage size
- Positive action locking mechanism
- Power washable
- Grease-free maintenance

Model 6250 Standard Features

- Conforms to BS EN 1865

Model 6251 Standard Features

- 4-inch (10 cm) front caster wheels

Model 6252 Standard Features

- Innovative Stair-TREAD system
- Extendable upper control handle
- 4-inch (10 cm) front caster wheels

Optional Features

- Two-piece molded ABS seat
- One-piece traditional vinyl seat
- Colour-coded patient restraints⁴ (plastic buckles)
- Nonabsorbent patient restraints⁴ (metal buckles)
- Polyester patient restraints⁴ (metal buckles)
- Patient head support (model 6252 only)
- O₂ bottle holder
- Vinyl head support
- Polypropylene head support
- IV clip
- Foot support

Specifications¹

Models	6250	6251	6252
Height	36 in. (91 cm)	36 in. (91 cm)	37.5 in. (95 cm)
Width	20.5 in. (52 cm)	20.5 in. (52 cm)	20.5 in. (52 cm)
Depth²	28 in. (71 cm)	28 in. (71 cm)	28 in. (71 cm)
Folded Depth	8 in. (20 cm)	8 in. (20 cm)	8 in. (20 cm)
Weight	20 lb. (9 kg)	23 lb. (10 kg)	31.5 lb. (14 kg)
Maximum Load³	500 lb. (226 kg)	500 lb. (226 kg)	500 lb. (226 kg)

Warranty

- One year parts, labour and travel or two years parts
- Lifetime on all welds*
- Seller's warranty hereunder shall not apply if: (i) Goods are not used in accordance with instructions or are used for a purpose not indicated on the labelling, (ii) any repairs, alterations or other work has been performed by Buyer or others on such item, other than work performed with Seller's authorization and according to its approved procedures; or (iii) the alleged defect is a result of abuse, misuse, improper maintenance, accident or the negligence of any party other than Seller. The warranty set forth herein is conditioned upon proper storage, installation, use and maintenance in accordance with applicable written recommendations of Seller. The warranty furnished hereunder does not extend to damage to Goods purchased hereunder resulting in whole or in part from the use of components, accessories, parts or supplies not furnished by Seller.

Seller's sole obligation shall be to repair or replace, at Seller's option, any defective component or item and pay transportation expenses for such replacement.

Extended warranties available.

*7-year service life.

¹Dimensions are measured from the outermost edges of the main frame. Specifications are rounded to the nearest whole number. Conversions are calculated before rounding.

²Depth dimensions are measured with extendable handles retracted.

³Maximum load capacity is total weight distributed in accordance to basic human anatomy. EMTs must consider the weight of the patient, equipment and accessories when determining the total weight on the product.

⁴Patient restraint set includes two cross-chest restraints and one ankle restraint. Cross-chest restraints can also be used in a chest and lap restraint configuration.

In-service video included with every order.

Stryker reserves the right to change specifications without notice.

This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Stryker**. All other trademarks are trademarks of their respective owners or holders.

The products listed above are CE marked according to the Medical Device Directive 93/42/EEC.



Literature number: MTX476EN Rev A
MTX-RRD/GS 01/11