


Rugged-X

restraints



Rugged-X restraints

Rugged-X restraints are Stryker's first dynamically crash tested coated webbing restraints, certified to meet BS EN-1789 crash safety standard.¹ With a durable, cleanable, polyurethane coating, Rugged-X restraints give you the power to put your focus on what matters most in an emergency - your patient, not soiled equipment.



Fluid resistant

Rugged-X restraints are equipped with a polyurethane coated webbing designed to support cleanability. The bright yellow coated webbing can help caregivers to see the restraints on the cot, identify soils, and expedite cleaning time.



Cost of ownership

Easy to clean material can help reduce the need for multiple backups and the need to discard hard-to-clean restraints.



Compatibility options

Rugged-X restraints are compatible with Power-PRO XT, Performance-PRO XT and Power-PRO TL ambulance cots.



Fluid resistant patch

Our engineered stitching/RF weld design provides strength while protecting stitching from hard to clean soils.

Reliable. Durable. Resistant.

Reliable.

Dynamically crash tested to meet BS EN-1789 crash safety standard¹.

Light color helps caregivers identify soils on the restraint webbing.

Durable.

Polyurethane coating protects underlying webbing from damage due to abrasion.

Rated for extreme temperatures between -41°C to 54 °C.

Resistant.


Soil resistant coated webbing makes Rugged-X easy to clean.

Cleanable polyurethane coating can help to expedite the time to clean and dry restraints.



Rugged-X


restraint part numbers



Full Set
6506-001-430



Shoulder (one each)
6506-001-401




Double buckle strap
6506-001-402



Buckle and Tongue Strap
6506-001-403



Thigh
6506-001-404



Ankle
6506-001-405

ProCare Services

Every day, you count on your medical equipment to perform at its best. With ProCare Services, our highly trained technicians help ensure your equipment is ready to perform when needed. Our goal is to make life easier for your staff and clinicians with service plans designed to help you get the most out of your investment. Choose the service plan that meets your needs.

ProCare Protect

A complete service solution designed to maintain your equipment with annual preventive maintenance while protecting your investment with an unlimited number of repairs, including parts, labour and associated travel.

ProCare Prevent

Supports your preventive maintenance initiatives with an annual inspection of your equipment, including parts labour and associated travel.

Flex Financial Program

Our financial programs provide a range of smart alternatives designed to fit your organisation's needs. We offer flexibility beyond a cash purchase with payment structures that can be customised to meet budgetary needs and help to build long-term financial stability. Contact your account manager for more information.

1. Meets BS EN-1789 dynamic crash standard with Power-LOAD cot fastener system used with Power-PRO XT ambulance cot, Power-LOAD cot fastener system used with Performance-PRO XT ambulance cot, Performance-LOAD cot fastener system used with Power-PRO XT ambulance cot, Power-PRO TL ambulance cot with fastener.

Medical EMS

This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Flex Financial, Performance-LOAD, Power-LOAD, Performance-PRO, Power-PRO, ProCare, Rugged, Stryker. All other trademarks are trademarks of their respective owners or holders.

The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

This material is not intended for distribution outside the EU and EFTA.



Mkt Lit-1473EN Rev A.6
DFS 07/2017
2017-14364

Copyright © 2017 Stryker
www.stryker.com

Stryker's Regional Headquarters Amsterdam
Herikerbergweg 110
1101 CM, Amsterdam
The Netherlands
stryker.com/emea